



BA Series Micro Balances

| Balance

☰ Details Info


Features :

Easy read Large backlit GRAPHICAL LCD display with A.E.P, (Advanced Eye Protection) (NEW)
 Standard RS-232 C interface I Ps2 output.
 Built-in Internal Calibration with function like E1 standard weight.
 Initial Calibration.
 Auto calibration according to change in ambient temperature.
 User selected time calibration.
 Various application other than just simple Weighing like
 %
 Filling
 Check weighing
 Statistic
 Improved Repeatability for better results.

Specifications

MODEL	BA 27	BA 56

CAPACITY	2 gm	5gm
READABILITY	0.1 µg	1 µg
REPEATABILITY(+/-)	0.4 µg	1 µg
LINEARITY(+/-)	0.5 µg	5 µg
PAN SIZE(mm/inch)	50-30 Ø	
RESPONSE TIME	8-10 sec.	
CALIBRATION	Perfect Auto Internal time & temp Controlled smart Calibration	
UNITS OF MEASURE	G,mg,ct,gm,GN,mo,oz,dwt etc.	
OPERATING TEMPERATURE	15°C to 35°C	
POWER SUPPLY	AC Adapter 230V or 115 V+/-20%50-60Hz	
WEIGHING CHAMBER	Ø 90 x 90 mm	
NET WT/GROSS WT	10.20 kgs / 14.70 kgs	
PACKING SIZE	565 X 565 x 355 / 22.2" x 22.2" x 14"	

 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
 Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved