



BD605 - Semi Micro Balance

| Balance

☰ Details Info

Semi Micro Balance

acculab Semi-Micro Balances offers exactly right performance and convenient features that you need to accomplish for the most sophisticated weighing tasks. you can weigh just about anything from really light weight samples to heavy weight objects and If you want to do more than "Basic Weighing" the 7 application Programs will let you easily perform the

Following :

- Piece Counting
- Percentage (%) weighing
- Check weighing
- Filling
- Statistics
- Recipe
- Density measurement for solid and liquid [with optional eqp & software]

Additional equipment

PRINTER CPR 02	ADDITIONAL LCD DISPLAY
PC KEYBOARD	WEIGHTRACK COMPUTER SOFTWARE
VIBROSORB-ANTI VIBRATION TABLE	CELIBRATION WEIGHT E2,F1,F2,M1(NABL,NPL,RRSL,STANDARD)

Specifications:

MODEL	BD605
CAPACITY	60/210gm
READABILITY	0.01/0.1 mg
REPEATABILITY(+/-)	0.02/0.1 mg
LINEARITY(+/-)	0.04/0.2 mg
PAN SIZE(mm/inch)	85Ø/3.3"Ø
RESPONSE TIME	3-5 sec.
DISPLAY	Backlit LCD Graphical Display
CALIBRATION	Perfect internal celibration
UNITS OF MEASURE	G,mg,ct,GN,mo,oz,dwt
TARE RANGE	Full
OPERATING TEMPERATURE	5°C to 40°C
HOUSING DIMENSION(WXDXH)	185X470x280mm/7.3"x18.5"x11"
WEIGH APPROX	9.3kg/20lb
POWER SUPPLY	AC Adapter 230V or 115 V+/-20%50-60Hz



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: