



BH Series (1g/0.1g; External Calibration)

| Balance

☰ Details Info

High Capacity Electronic Balance

Features:

- 1. High precision load cell.
- 2. Rugged Plastic Housing.
- 3. Big size Stainless Steel Platter.
- 4. High contrast easy to read LED display.
- 5. Built-in rechargeable battery last up to 7 days.
- 6. Height adjustable feet.
- 7. Full capacity subtraction
- 8. Stabilization time 1.5 seconds
- 9. Panel simple Parana facile operation
- 10. Overload Protection
- 11. Below Balance Weighing Facility
- 12. Checking Weighing, piece counting function.
- 13. RS232 Interface.

Model	BH60001	BH100001	BH150001	BH200001	BH300001	BH15000	BH30000
Capacity	0-6000g	0-10000g	0-15000g	0-20000g	0-30000g	0-15000g	0-30000g
Readability	0.1g	0.1g	0.1g	0.1g	0.1g	1g	1g
Repeatability	±0.2g	±0.2g	±0.2g	±0.2g	±0.2g	±2g	±2g
Linearity	±0.2g	±0.2g	±0.2g	±0.2g	±0.2g	±2g	±2g
Pan size	335*235mm						
Dimensions	370*370*125mm						

Net Weight	5000g
Power Supply	DC 7.5V/350mA /Rechargeable Battery 4V/4AH



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved