



## BP-D series ( 0.1g/0.01g; External Calibration) with bottom loading

| Balance

☰ Details Info

BP-D series density determination electronic balance( 0.1g/0.01g; External Calibration) with bottom loading

Features:

1. High precision load cell.
2. Rugged Plastic Housing.
3. Big size Stainless Steel Platter.
4. High contrast easy to read LED display.
5. Built-in rechargeable battery last up to 7 days.
6. Height adjustable feet.
7. Full capacity subtraction.
8. Stabilization time 1.5 seconds.
9. Overload Protection.
10. Below Balance Weighing Facility.
11. Checking Weighing, piece counting function.
12. Can be contain bubble level.

Specifications:

Model	BP50001D	BP50002D
Capacity	0-5000g	0-5000g
Readability	0.1g	0.01g
Repeatability	±0.2g	±0.02g

Linearity	±0.2g	±0.02g
Pan Size	143*192mm	143*192mm
Container	Ø20*24cm	Ø20*24cm
Basket Size	Ø16*27cm	Ø16*27cm
Dimensions	390mm x 250mm x480mm	390mm x 250mm x480mm
Net Weight	4800g	4800g
Power Supply	110V/220-240V AC	110V/220-240V AC



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa  
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

[sales@acculabusa.com](mailto:sales@acculabusa.com) [middleeast@acculabusa.com](mailto:middleeast@acculabusa.com)



Phone Support:

---

© Copyright 2014 by ACCULAB\_USA - All rights reserved