



IRS Infrared Sterilizer

| Biosafety Accessories

☰ Details Info

Features

A mini-type sterilizer, it adopts inside radiation components and numeric display temperature control technique. It is made of stainless steel , its characteristic includes low weight, small dimension , quick sterilization and safety, energy-saving , clean.

Specification:

Input voltage	220 V
Sterilization temperature	280°C±10°C
Frequency	50HZ
Average consumption	200w
Protection type	1
External dimension	170×100×170 mm
Net weight	1.85kg

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

✉
Email Us At: