



MVM-50 Multi-Tube

| Vortex

☰ Details Info

Specifications

1. Speed range: 500-2500rpm
2. Working diameter: 3.6mm
3. 50 samples

MVM-50 is mix one can handle up to 50 samples. Speed is adjustable from 500 to 2500 rpm with electronic speed control. Built-in timer is adjustable from one minute to 99h59min operation. Various tube sizes can be mixed by simply changing the foam tube racks.

Pulsing Mode

Programmable pulsing mode allows user to adjust the pulse-on and pulse-off times between 1~59 seconds in 1 second intervals. This feature enhances the vortex action by creating a more vigorous mix.

Features

1. Can handle 50 samples
2. Two modes: continuous type of work or Programmable pulsing mode
3. Microprocessor control
4. LCD display, and easy to use the panel
5. The variable speed microprocessor control provides consistent uniform mixing action

6. Independent LCD displays for speed and time allow operator to view both settings at once

Technical parameters:

Speed range	500-2500rpm
Speed accuracy	± 50rpm (at 2500rpm detection)
Timing range:	1min-99h59min
Mixing Orbit	3.6mm
Maximum load	4.5kg
Dimensions	426x250x480 mm
Weight	21.5kg
Power	230V or 120V, 50/60Hz, 60W

Optional Accessories:

Type	Detail	Dimension (mm)
MVM-B01	50 x Φ10mm tube frame	245×132×45 mm
MVM -B02	50 x Φ12mm tube frame	245×132×45 mm
MVM -B03	50 x Φ13mm tube frame	245×132×45 mm
MVM -B04	50 x Φ16mm tube frame(15ml tube)	245×132×45 mm
MVM -B05	15 x Φ25mm tube frame	245×132×45 mm
MVM B06	15 x Φ29mm tube frame(50ml tube)	245×132×45 mm
MVM -B00	Tray Pad Set (Upper & Lower)	305×178.5×25 mm

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved