



UPW-L Series

| Water Purification

☰ Details Info

Water purifier (Ultra-pure grade)

Features:

Professional Water purifier produces ultrapure water.
 Widely used in hospital, research lab, industry, university etc.
 Automatically power on and power off according to the water Inflow or water outflow.

LCD Display	Power Process indication State indication
	Power-off protection
Self-Safety	Auto water-level control
	Low Input Pressure protection
	Auto restart Visual/Audio alarm

Specification

Model	UPW-30L	UPW-80L	UPW-120L	UPW-200L
Capacity	30L	80L	120L	200L
Ultra pure water				
Output water quality	≥15MΩ.cm (25°C Ambient)			
TOC(RO water)	1-10 ppb			
Bacteria	<1 CFU/ml			
PH	7			
Index				
Dimension (mm)	450*250*650	420*360*920	500*400*1150	700*600*1600
Weight	20 Kg	40 Kg	50 Kg	70 Kg
LCD Display	Power Process indication State indication			
	Power-off protection			
Self-Safety	Auto water-level control			
	Low Input Pressure protection			
	Auto stop	Visual/audio alarm	Visual/audio alarm	Visual/audio alarm

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved