



## ARJ-3500

| Anaerobic jar

☰ Details Info

### Product Description:

Transparent Anaerobic Culture Jar of 3.5.lit. capacity which is made of Polycarbonate.

The lid of the ARJ is made of Aluminum pressure die Casting fitted with Pressure Valve,

Vacuum Valve, Safety Valve & compound Pressure cum vacuum gauge, three finger aluminum clamps and sealing "O" ring.

This is recommended only where inside air of ARJ is removed thru vacuum pump and CO<sub>2</sub> is filled in the jar thru CO<sub>2</sub> gas cylinder to generate anaerobic condition in the ARJ.

Also supplied with plain aluminum pressure die casted lid which can be fitted when readymade gas packs are used to generate the anaerobic conditions inside the ARJ.

When readymade gas packs are used, there is no need of lid fitted with dial gauge etc., for extra life of equipment.



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAoffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa  
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: