



Armrest

| Biosafety Accessories

☰ Details Info

Small dimensions and easy to handle, with a adjustable height of 20 mm . Its design conforms to the airflow grill of biosafety cabinet, which avoids influencing the airflow pattern during operation . It will help reduce the fatigue of operators' hands.

Size:120 mm×80mm×82mm



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAoffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved