



BP series 0.01g

| Balance

☰ Details Info

Features:

1. High precision load cell, compact design.
2. Full plastic housing.
3. Stainless steel Plstter.
4. Super bright LCD display with backlight.
5. Mains adapter supplied as standard.
6. Height adjustable feet.
7. Full Capacity subtraction.
8. Simple 5 key panel for easy operation.
9. Stabilization time 1.5 seconds.
10. Overload Protection.
11. Below Balance Weighing Facility.
12. Selectable measure units, g,oz,ct.
13. Checking Weighing, Piece counting function.

Specification;

Model	BP30002	BP40002	BP50002
Capacity	0-3000g	0-4000g	0-5000g
Readability	0.01g		
Repeatability	±0.02g		

Linearity	±0.02g
Platter Size	155*143mm
Dimensions	280mm x 180mm x80mm
Net Weight	1200g
Power Supply	DC 6V-9V(AAA battery *4pcs)



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved