



Blood Centrifuge Balance

| Balance

☰ Details Info

Blood Bank Centrifuge Counting Balance

It is a specially designed electronic differential scale for a particular use in blood bank. It measures two weights simultaneously and displays weight difference instantly. The first step in plasma separation from the whole blood is loading of blood bags in high speed centrifuge machine. Before loading the blood bag, the weights of the two loads (each blood bag + bucket) have to be balanced in such a way that their weights should be the same or very close to each other. This is a vital step of plasma separation process to avoid damage to centrifuge machine. By conventional methods the process of balancing is time consuming, inaccurate and cumbersome. This equipment makes the whole process of balancing execute shortly and accurately. It is the only way to save expensive centrifuge machine from damage

Features

- Accurate, simple and easy to use
- With two independent weight sensors, which displays individual weight of each bucket with high degree of accuracy
- Can weigh up to 5000grams
- Indicates the weight of each pan separately
- With a bright LCD display
- Displays the difference of the two weights
- The bright LCD glows when both the weights on the pans are balanced(Equal)
- Audio overweight(>5000gms.) alarm

Hassle free operation by built-in battery backup.
 Leveling screws for leveling the equipment
 Water bubble level indicator
 External calibration screws for quick and efficient calibration
 Every equipment accompanies a calibration certificate of a standardized lab
 Built-in replaceable fuse for extra safety of the equipment
 Power 220-240V AC Supply

Specification:

Weight Measurement Range	0-5000 grams
Weight Difference Range	5000 grams (max)
Weight Sensors	Load Cells
Weight Display Panel	4 digit LED display
Equal weight Indi.	RED BRIGHT LED
Size of single pan	206 X 138mm
Accuracy	+/- 1 gram
Voltage Requirements	220-240 Volt A.C.
Battery Backup	Yes
Weight	4.5 kg. approx.
Dimension	310x280x100mm



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
 Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: