



CC-31 Le.Chatelier Cement Steam Bath

| Mining Industry

☰ Details Info

Specifications:

- 1.Le. chatelier cement steam bath
- 2.test soundness of cement
- 3.operation easy
- 4.best steam tank supplier
- 5.3 years warranty

Equipment summary:

For curing the specimen for the test of soundness of cement with both Le. Chatelier and pat methods. Two heaters are fitted to maintain an accurate temperature control in the water.

Technical parameter:

~

Model	LCC-31 Le. chatelier Cement Steam Bath
Effective volume	31L
Heating time	3h
Heaters power	4kW
Net weight	≈20kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved