



Formalin Fumigation Sterilizer

| Biosafety Accessories

☰ Details Info

Specification:

Model	FORS	AHCN
Voltage/Frequency	Voltage/Frequency: 220V/50 Hz	
Powe	180w	
Volume	100 ml	
	It used for sterilization of the safety cabinet before changing filte	It used for toxicity elimination of the Formalin

📍 Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

✉ Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

☎ Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved