



HSG10-50- Grinder

| Grinder

☰ Details Info

Product Description:

Summary of Use: It applies in grinding all kinds of food, soil, medicine, ore and mineral substance in the fields of industry, agriculture, industrial mining, medicine and health, colleges and universities, coal and geological etc.

Characteristics:

- 1.The grinding room and cutter blade are made of imported stainless steel for the accuracy of analyzing treatment of experimental article.
- 2.It has advantages of small cubage, high grinding efficiency, easy operation, beautiful shape and wide use.

Major technical parameter:

Working voltage (V)	$\pm 10\%$ 220V... 50Hz
Power (W)	180

Speed (r/min)	10000
Size of grinding room (mm) (diameter×high)	Ñ,80x40
Crushing capacity (≤g)	50
Fineness of grinding (mesh)	40-200
Overall dimension (mm) (diameter×high)	Ñ,110x210
NW (kg)	2.5
GW (kg)	2.7



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved