



KVT-1007

| Kinematic Viscometer

☰ Details Info

Specifications:

Power supply	AC 220V±10%, 50Hz±5%
Heating power	two grades, 1000 W+650 W
Stirring motor	6 W, 1200 RPM
Temperature range	ambient temperature ~100°C
Temperature controlling accuracy	±0.1°C
Constant temperature bath	20 L; two layers water bath;

Operation environment	(1) Ambient temperature: -10°C~+35°C (2) Relative humidity: <85%
Temperature sensor	Industrial Platinum Resistance; Pt100
Capillary Viscometer: 7pieces in a group. Inner diameters of them are 0.6, 0.8,1.0, 1.2, 1.5, 2.0, 2.5mm respectively (inner diameters of 5.0 and 6.0 mm areoptional);wer consumption: not more than 1800 W	



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved