



KVT-2004 - Automatic

| Kinematic Viscometer

☰ Details Info

Features:

Microprocessor controls the total test, with menu for operation.
 Automatic control the water bath, timing, viscosity calculation, capillaryviscometer washing and drying, print the report.
 ·Automatic self diagnostic, remind where, how in faults, so that to do it efficiently.

Specifications:

Power supply	220V±10%, 50Hz
Capillary viscometer	U-viscometer
Temp. accuracy	±0.1°C
Timing accuracy	±0.1s
Max time	999.9s

Max consumption power	2000W
-----------------------	-------

Note: standard configure is with 4 capillary viscometers, ID 0.8mm, 1.0mm, 1.2mm and 1.5mm, if you want others capillary viscometer, please order it additional.



Head Office : 1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USA Office: 128 W SUFFLOK AVE, CENTRAL ISLIP, 11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved