



L JC- Series

| Mining Industry

☰ Details Info

L JC- Series Coal sample preparation lab jaw crusher

Features

1. Jaw plates with BTMcr20 material are reversible for extended life
2. Patented design for easy and quick jaw gap adjustment provides good control of product sizing while minimizing downtime and labor service requirements.
3. Hermetic structure improving the operation environment and reducing sample contamination.
4. Improved safety with protection shielding and emergency button.
5. Access to the chamber through the front door makes it convenient to clean.

MODEL	L JC-500 Jaw Crusher	L JC-1000 Jaw Crusher
Size of Feed Inlet	100×60 mm	150×125mm
Feed Size	≤50mm	≤100mm
Discharge Size	≤3~13mm adjustable	≤6~38mm adjustable
Moisture Adaptability	≤12%	≤12%
Approx. Hardness of Feeding Sample	HRC48-55	HRC48-55
Throughput	200~500kg/h	500~2000kg/h
Power Supply	3 phases, 380V	3 phases,380V
Power of Electric Motor	1.5kW	4kW
Net Weight	201kg	378kg
Dimensions(L×W×H)	875mm×460mm×765mm	1046mm×600mm×1090mm



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved