



LVS-1300 Viscometer

| Tile Industry

☰ Details Info

**Specification:**  
**13,000,000 cp Engler viscometer**

<b>Measurement Range</b>	1-13×10 <sup>6</sup> mpa.s
<b>Display</b>	LCD
<b>Spindles Specification</b>	No. 1, 2, 3, 4, 5, 6, 7( No. 0 is optional) No. 0 test 1-100 mpa.s Note: The bigger the number of the Spindles, the bigger the measuring value
<b>Rotating Speed of Spindles</b>	0.3 / 0.6 / 0.5 / 3 / 6 / 12 / 30 / 60 rpm/min. with 1, 2, 3, 4# Spindles 0.5 / 1 / 2 / 2.5 / 4 / 5 / 10 / 20 / 50 / 100 rpm/min. with 5, 6, 7# Spindles
<b>Interface</b>	Printer Interface PC Interface Temperature Sensor Interface

<b>Automatic Switch</b>	free selection of proper rotating speed or Spindles number.
<b>Stable Cursor</b>	Strip-type cursor fully charged displayed when reading basically stable
<b>Measurement Tolerance</b>	± 1%
<b>Power</b>	AC 220V ± 10% 50 Hz ± 10%
<b>Working Condition</b>	Temperature 5-35 degree Relative humidity no more than 80%
<b>Attachments</b>	No. 1, 2, 3, 4, 5, 6, 7 Spindless (No. 0 Spindles is optional, if need, please contact with me.) Adapter Protection Frame Foundation Plate Up-Down Frame English Manual Hexagonal Spanner Solid Wrench
<b>Outline Dimension(mm)</b>	370×325×280 mm
<b>Packing (mm)</b>	400×360×310 mm
<b>Net Weight (kg)</b>	6.8 kg
<b>Gross Weight (kg)</b>	9 kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N  
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: