



LWD1 Lab Waste Disposer

| Biosafety Accessories

☰ Details Info

Features:

Lab Waste Disposer is a dry sterilizing .Lab Waste Dispose with biosafety guarding function. The system is used for disposing the on-site solid wastes caused by a little high-risk infectious or pathologic inspection in hospitals or labs, including various infectious medical wastes, especially including the excision body of human and animals; microbiologic incubation, PCR experiment wastes; needle and syringe of injection, etc. The wastes can be treated as the common wastes after disposed by this system.

Main characteristics:

1. Interior structure adopts aviation high-temperature material to avoid the spare parts appearance aging caused by high temperature during the sterilization. Exterior is made of fully stainless steel, easy to clean and sterilize.
2. There is the high efficient active carbon filter equipped at the exhaust vent of the cooling pipe, which can prevent the leakage of aerosol and particular smell.
3. Adopts interlocking control to avoid the door random opened before the sterilization is finished, also to prevent the diffusion of aerosol and germs.
4. It has inner high efficient thermic insulant to assure the surface temperature is no more than 45°C during sterilization and heating process and avoid human hurt caused by wrong operation.
5. The little volume waste-liquid is collected in a specified bottle to assure the safety of environment and personal.
6. During sterilizing, microcontroller controls the temperature and time to assure the sterilization effects.

Specification:

| | |
|----------------------|---|
| Dimensions(mm) | 360*260*230 |
| Process Capacity(L) | 0.6 L |
| Grounding resistance | ≤ 0.1Ω |
| Power supply | 110~240V/50HZ/60HZ |
| Process temperature | above 170 centigrade |
| Process time | above 170 centigrade, maintain over 1.5 hours |

📍
Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved