



## MID-2000

Industrial Microscopes

## **■** Details Info

## Features:

MID-2000 DIC Industrial Microscope has been experimented and advanced many times by our professional staff, has steady performance and good effect, is ratified by several top 500 companies, has reached the domestic leading standard. Welcome distributors to consult.

Specification		
Optical System	Infinity Optical System F=200mm	•
Viewing Head	Siedentopf trinocular head, inclined at 30°, interpupilary 50~75mm	
Eyepiece	WF10X/20mm	
	WF15X/12mm	

710		Accurabusaivi		
Objective	Focus Length 45mm	Plan Achromatic Objective 4X/0.10 WD=17.41mm	•	
		Plan Achromatic Objective 10X/0.25 WD=16.29mm	•	
		Plan Achromatic Objective 20X/0.45 WD=1.06mm	•	
		Plan Achromatic Objective 40X/0.65 WD=0.70mm	0	
		Plan Achromatic Objective 50X/0.75 WD=0.562mm	0	
		Plan Achromatic Objective 60X/0.75 WD=0.579mm	0	
	Focus Length 95mm	M Plan APO HL Objective 2X/0.05 WD=20.00mm	0	
		M Plan APO HL Objective 5X/0.12 WD=22.00mm	0	
		M Plan APO HL Objective 10X/0.30 WD=34.00mm	0	
		M Plan APO HL Objective 20X/0.30 WD=33.00mm	0	
		M Plan APO HL Objective 50X/0.5 WD=20.70mm	0	
Negopiose	Triple	•		
Nosepiece	Quintuple nosepiece with the DIC Jack			
Stage	Double layers stage 350×255mm, moving range: 200×200mm			
Filter	Green, Blue, Neutral			
Focusing	Coaxial coarse and fine with rack and pinion, fine scale value 0.002mm			