



UAZ-20VET

| Urine Analyzer

☰ Details Info

Features:

1. Point – of –care urinalysis with palm –sized feather weight units
2. Bluetooth technology for wireless printing and data transfer
3. LCD colored touch –screen
4. Single button easy operation
5. Recharge , transfer data and upgrade ,all with Micro USB cable
6. Compatible with URIT 10G,11G AND 14 G urine reagent strips

Specifications:

Model	UAZ-20VET
Parameters	PH, Nitrite, Specific Gravity, Blood , Urobilinogen , Ketone, Micro albumin , CalciumBilirubin, Leucocyte, Protein,Glucose, CreatinineAscorbic Acid , Micro albumin , Calcium
Measurement Principle	Dual wavelength reflectance method

Measurement Wavelength	4 LED wavelength (470,550,620,720 nm)
Sample Preparation	Manual dipping
Reaction time	60 seconds per test
Data memory	Up to 1000 samples
Display	LCD Colored touch – screen
Strips	URIT 10G/11G(14G optional)urine reagent strips
interface	Micro USB, Bluetooth
Work condition	80%~C, RH:20%°Environmental temperature 5-40
Power supply	Power adapter 9V DC,1.5A
Dimension	120(L)X83(W)X31(H) mm (analyzer) 110(L)X80(W)X38(H) mm (printer)
Weight	180g(analyzer), 180g (printer)



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved