



UAZ-50VET

| Urine Analyzer

☰ Details Info

Features:

1. Cold light source facility with long life and high stability
2. Internal thermal printer
3. Easy operation, cost-efficiency and high quality
4. Rs232 serial interface to host PC/LIS, Barcode reader optional
5. ·URIT UQ-10 & UQ-11 Urinalysis Control
6. URIT 10G、11G
7. urine reagent strip
8. Complete urine analysis with Ascorbic Acid

Model	UAZ-50VET
Parameters	pH, Nitrite, Specific Gravity, Ascorbic Acid, Blood, Glucose, Bilirubin, Urobilinogen Ketone, Leucocytes, Protein,Calcium
Measurement Principle	Dual wavelength reflectance method

Measurement Wavelength	3 LED Wavelength (550,620,720nm)
Sample Preparation	Manual dipping
Reaction time	60 seconds per test
Throughput	Normal mode, 60 samples /hour Continuous mode, 125 samples /hour
Memory capacity	5000 results
Display	LCD Colored touch – screen
Strips	URIT 10G /11G urine reagent strips
interface	,(Serial port(RS232), parallel port (For external printer PS/2(For barcode reader
Work condition	, RH: 20%~80%°C Environmental Temperature: 15-30
Power supply	Adopter 12V DC, 4.0A
Dimension	290mm x 200mm x 100mm
Weight	2kg



Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N
Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE



Email Us At:

sales@acculabusa.com middleeast@acculabusa.com



Phone Support: