



WM-1

 | Whiteness Meter

[Details Info](#)
Specification :

Measuring range	0 - 199
Whiteness formula	blue whiteness WB = R457
Light source	LED light source
Display	LCD liquid crystal display
Measuring bore diameter	φ 15
Minimum Indication	0.1
Measurement repeatability	0.2
Indication error	≤ ± 1

Test conditions: compliance with the provisions of GB/T3978 45/0

Accessory :

Silver crucible 30ml	content: 99.99% capacity: 30 ml
Platinum crucible 30ml	content: 99.99% capacity: 30 ml

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

Phone Support: