



XYS100M

| X-Y Positioning stage

☰ Details Info

Precision Manual Stage X-Y

Specifications

Precision Manual Stage for Positioning for Semiconductor machine, Telescopic, Medical equipment, and many other positioning equipmen

Stage size 100 mmX100mm

Travel 70 mm in Xand Y position

Travel per knob rotation 20 mm

Load capacity 5 kg

include all bolts and other fixtures required to fix the stage into a workstation base

Stage constructed by durable material that will not rust in case splashed with water nor corroded from common laboratory chemical such as acetone and alcohol

📍
Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa Office
: warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

✉
Email Us At:

sales@acculabusa.com middleeast@acculabusa.com

☎
Phone Support:

© Copyright 2014 by ACCULAB_USA - All rights reserved